

The thermal conductivity analyser **HLC A310** for the production-area and stock-receipt measures thermal conductivity and thermal resistance of 300*300mm platform insulation materials.

Like all **HLC**-devices this stand-alone unit has to be connected on a standard 230Vac-socket only and needs no water-connection and stands out for easy handling also.



According standards

EN12667:2001
EN1946-3:1999
ISO8301

the instrument measures the thermal transfer properties under steady-state conditions of
| 300*300mm platform insulation materials
| with a thickness of 15...100mm

in meas-ranges of ...


| thermal resistance R up to 20m²*K/W
| transfer factor τ resp. thermal conductivity λ 5...200mW/(m*K)

By electromotoric positioning

- not only soft insulation materials with automatic distance-positioning
- but also compact foam materials by using the pressure-position-mode on 225N

| can be analyzed under exact repeatability conditions.

The relative meas method was realized with two sensitive heat flow sensors (symmetrical configuration) for a mean temperature of 23°C. With a non metering edge-area of each 75mm up to 100 mm thick insulation materials can be measured at ambient temperature 23°C.



Symmetrical configuration

H'' heat plate U'' cool plate
 H', H'' heat flow meter
 S sample
 T_m controlled mean temperature

Construction

The heat and cool copper plate is exactly temperatured with a PI-regulated peltier-element on 15°C resp. 31°C. The peltier-units of both secondary sides are against-temperatured with a powerful liquid circulation.

On both plate surfaces are sensitive heat flow sensors adapted..

Through a servo-electronic the lower cool-plate with inserted sample is driven against upper heat-plate. Fixed in the center of the springy cool-plate a linear-measurement-system for automatic thickness registration is installed, a pressure switch under cool plate limits the servo-controlled pressure at 225N.

An independent **HLC A310**-electronic with power supplies for the temperatures control systems and sensor signal read-out is integrated. On a connected PC or Laptop *WINDOWS*-software **WinHLT#** calculates in Online-mode all meas values, which will be displayed in tables and graphic diagrams on desktop.

Features


With the Lambda-range until 200mW/(m*K) all conventional insulation materials are measurable, also plates with metallic cover. In consequence of the thermal hemispherical emittance > 0,8 this analyser is for low density materials especially qualified.


Calibration of the relative measurement is factory-set with several calibration-samples determined per **ISO8302**. With own calibration material, the user can always rapidly and simply test the adjustment.

The instrument needs no special requirements and can be used continuously. When a several day stop occurs, a Standby-mode with switched-off power-consumption can be selected.

WINDOWS-Software WinHLT#

With **HLC A310** supplied WINDOWS-software **WinHLT#** runs under **MS-WINDOWS XP / WINDOWS7**. It records, calculates and archives meas data.

Testreport according EN12667:2001 section 9	
a) Test determination:	Determination with heat flow meter methods
b) Kind of apparatus:	Single-specimen symmetrical configuration
c) Identification of apparatus:	HESTO-Lambda Control A206 - SN: 619
d) Specimen standard:	HESTO
e) Specimen no.:	0
f) Name and partment identification of specimen:	XP82030.8
g) Description of specimen:	
h) Conformity of specimen:	
i) Preparation according standard:	
j) Specimen thickness:	0,0505 m - automatic measured
k) Preparation according standard:	
l) Density of prepared specimen:	
m) Rel change of mass during pre-treatment:	
n) Rel change of mass during measurement:	
o) Middle temperature difference between specimen surface:	16,16 K
p) Middle-Temp.:	23,00 °C
q) Density:	(15,84 + 16,88) / 2 = 15,86 W/m²
r) Thermal Resistance:	1,92 m²/KW
s) Thermal Conductivity:	0,0466 W/(m*K) [λ] = 9999,990 mW/(m*K)]
t) Date of test end:	12.16 / 02.09.2013
u) Duration of test:	230 min
v) Date of last sensor calibration:	-no calibration data-
w) Calibration-Spec.:	Thermal Resistance
x) Date of certificate:	
y) Certificate Inset.:	
z) cert. number:	
aa) Orientation of the apparatus:	horizontal, warmside up
ab) Vibration-resistant cover:	
ac) Meas graphic:	
ad) Meas deviation:	2,8 % = $\sqrt{2 \cdot 5^2 + 2 \cdot 5^2}$
ae) Conformity:	This test conforms to EN12667:2001 "Thermal performance of building materials and products - Determination of thermal resistance by means of guarded hot plate and heat flow meter methods - Products of high and medium thermal resistance"
af) Name of operator:	HESTO



With end of measurement the program automatically prints the test report, which is an important request of **EN12667:2001**; additional test declarations for this report will be insert in remanent text-fields. **WinHLT#** stores all these data as a data-set in a table-file on the hard-disc also.

The automatic run will be moderate with meas-menus and grafic diagrams, variable keys leading errorfree to several menus for different modes and parameters. All modes keep stored after meantime analyser or PC power-off also (means an interrupted measurement in this case will continue automatically).

So with unchanged meas-mode normal operation is restricted with inserting sample material into the measurement chamber and pushing the start-key.

WinHLT# includes several language-variants: actual german, english, french, spanish, italian and polish.

Prinzip and method of operation

After the specimen was insert and "start"-button was pushed the measurement chamber will be closed by electromotoric positioning while a linear-measurement-system (solution 0,05mm) automatically registrade the thickness. On the upper specimen surface a copper-plate is exactly temperatured on 31,0°C with a PI-regulated Peltier-element. An equal unit temperatures the lower specimen surface on exact 15,0°C. Two sensitive heat flow sensors are integrated on both plate surface.


HLC A310 transfers the meas-values to the PC. With the adjust-parameters (stored on hard-disc) **WinHLT#** calculates as follows:

Sign	Size	Unit
R	thermal resistance	m²*K/W
λ	thermal conductivity	W/(m*K)
A	sensitive area	m²
d	average thickness	m
T₁	temperature hot side	K
T₂	temperature cold side	K
φ	heat flow rate	W

$$R = \frac{A * (T_1 - T_2)}{\phi}$$

$$\lambda = \frac{\phi * d}{A * (T_1 - T_2)}$$

The necessary physical thermal balance (means steady-state condition) is ready, if both heat flux sensors deliver stable and similar signals. This status is indicated on the desctop and with a blinking LED "measurement" on **HLC A310** front. If operator finish measurement via select in menu now, chamber opens electromotorically and moves out specimen. **WinHLT#** print out the test report (as pdf and/or on a connected printer) and stores final meas data as a data-set on hard-disc.



3 year guarantee and manufatorer direct-service

The instrument needs no special requirements and is designed for continuous operation time. Produced in Germany near Frankfurt/Main, the reliability and durability is protected with a **3-year guarantee**. In service-case the manufacturer helps fast and direct.

Inspection and Cleaning Service

Quality assurance according to **EN ISO9000 ff.** is standard in the insulation industry. The accordance standard prescribes as an elemental condition a cyclic check with certification of the test equipments. Accomplish to this demand we offer for **HLC A310** a low-price and fast inspection/cleaning service inclusive acceptance-protocol with meas date, results and links to used reference materials.

Optional automatic Sample-Handling-System RSH3

With the automatism specimen-handler option **RSH3** the analyser **HLC A310** becomes a fully automatic meas-system. An industrial robot picks measured plates out of meas-chamber and puts it on a storage-stack, thereafter he gets next to be measured plate from a pick-up-stack and brings it into analyser. So one after another could be measured without any manual intervention stacked specimens between 15...100mm thickness in size 300*300mm. Specimens with Barcode on an adhesive label will be identified with a included reader-unit and last meas-data together with barcode-number and -text will automatically stored on end of measurement as a data-set in a file on PC harddisc.

General Specifications HLC A310

Construction accord. ISO8301	Single-specimen symmetrical configuration, heat flow meter on heat- and cold-plate each
Specimen thickness automat. Meassystem	15 to 60mm (according to EN1946-3:1999) incremental Linear-Measurement; Display Solution 0,1mm
Specimen size, weight	300*300 mm, max. 1,75kg
Case sensitive area	each 150mm
Width of non metering area	each 75mm
Measurement-Range	λ (T) 5...200mW/(m*K); additional calculation + display of λ_{10} value according EN10456 R up to 20m ² *K/W
Accuracy	< \pm 3% at 23°C ambient temperature
Repeatability	< \pm 1%
Mean temperature	23°C , \pm 1°C
Temperature hot plate	31,0°C < \pm 0,1°C
Temperature cool plate	15,0°C < \pm 0,1°C
Ambient temperature	23°C, \pm 1°C
Positioning	electromotoric, distance selectable from 15,0 to 100,0mm; selectable pressure positioning with 225N
Meas time	about 20 minutes for 20mm thickness;for thicker material up to several hours until steady-state condition
Warmup time	about 120 minutes after power-on
PC-Interface	USB2.0
Power Requirements	230V/50Hz about. 300W
Instrument Size	600* 750*600mm (B*H*D)
Noise	about 40dB
Instrument Weight	about 100kg
Delivery scope	analyser HLC A310 with power- and pc-interface-cable, a XPS/EPS test sample for cyclic check of analyser ("Internal control"), pc-software WinHLT# (pc is not part of standard delivery)